WOUND CARE

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1. "A Closer Look At Ultrasonic Debridement" Issue Number: Volume 23 - Issue 8 - August 2010 (/issue/2235), Author(s); Martin E. Wendelken, D.P.M., RN, Lee Markowitz, D.P.M., and Oscar M. Alvarez, PhD"
2. "Split Thickness Skin Graft Incorporation after Ultrasonic Debridement for Wound Bed Preparation in Diabetic Foot Wounds." Sean M. Betesh, D.P.M., Michael I. Gazes, DPM, MPH, Peter Blume, DPM, FACFAS.
3. "Ultrasonic Tangential Burn Excision Reduces Blood Loss." Abraham P. Houg M.D., Sylvia J. Petrone M.D., Robin A. Lee M.D., Christina Lee M.D., Michael A. Marano M.D.
4. "Surgical Debridement with SonicOne; an initial look at the bacterial load pre and post low frequency direct contact ultrasound." Daria Abolghasemi, D.O., Karen D. Szymanski, D.O., Michael Baruch M.D., Jamshed Zuberi M.D.
5. "Ultrasonic Debridement An Effective and Superior Method of Wound Management." John M. Hiebert M.D. Abstract presented at SAWC 2013.
6. "Ultrasonic Tangential Burn Excision Reduces Blood Loss." Abraham P. Houg M.D., Sylvia J. Petrone M.D., Robin A. Lee M.D., Christina Lee M.D., Michael A. Marano M.D.
7. "Applications of the Ultrasonic Bone Cutter in Spinal Surgery - Our Preliminary Experience." Al-Mahfoudh R, Qattan E, Ellenbogen JR, Wilby M, Barrett C, Pigott T. Br J Neurosurg. 2013, Jul 10.
8. "Can Ultrasound Debridement Facilitate Biofilm Removal from Diabetic Foot Ulcers?" Melinda Bowlby, D.P.M., Peter Blume, D.P.M. Podiatry Today. Volume 27 - Issue 8 - August 2014.
9. "Surgical Debridement with SonicOne; an initial look at the bacterial load pre and post low frequency direct contact ultrasound." Daria Abolghasemi, D.O., Karen D. Szymanski, D.O., Michael Baruch M.D., Jamshed Zuberi M.D.
10. "An overview of factors maximizing successful split-thickness skin grafting in diabetic wounds." Ryan J. Donegan, D.P.M., Brian M. Schmidt, D.P.M., Peter A. Blume. D.P.M. Diabetic Foot & Ankle.

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



## STUDI CLINICI

*Surgical Debridement with SonicOne® O.R.: an Initial Look at the Bacterial Load Pre- and Post- Low Frequency Direct Contact Ultrasound*

**Daria Abolghasemi, D.O.1, Karen D. Szymanski, D.O.2, Michael Baruch M.D.3, Jamshed Zuberi M.D.4**

### Affiliations

1. Resident Department of General Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey
  2. Fellow Section of Plastic and Reconstructive Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey, USA
  3. Chief Section of Plastic and Reconstructive Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey, USA
  4. Trauma and Critical Care Surgeon, Department of Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey, USA
- Abstract Presented at The Diabetic Limb Salvage Meeting, March 2016

*Split Thickness Skin Graft Incorporation after Ultrasonic Debridement for Wound Bed Preparation in Diabetic Foot Wounds*  
**Sean M. Betesh, DPM1, Michael I. Gazes, DPM, MPH1, Peter Blume, DPM, FACFAS2**

1. Resident, Department of Podiatric Surgery, Yale New Haven Hospital, New Haven, CT, USA
2. Assistant Clinical Professor of Surgery, Anesthesia, and Orthopedics and Rehabilitation, Yale School of Medicine, New Haven, CT, USA

Abstract Presented at The Diabetic Limb Salvage Meeting March 2016

*The Clinical Implications of a New Wound Debridement Device That Combines Low Frequency Direct Contact Ultrasound and Vacuum Aspiration*

**Mark S. Granick M.D.1, Michael Baruch M.D.2, Wayne J. Caputo M.D.3, Paul M. Glat M.D.4**

### Affiliations

1. Division of Plastic Surgery, Rutgers New Jersey Medical School, Newark, New Jersey, USA
2. Division of Plastic and Reconstructive Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey, USA
3. Clara Maass Medical Center, Barnabas Health, Belleville, New Jersey, USA
4. Division of Plastic Surgery, Division of Burn Unit, St. Christopher's Hospital for Children, Philadelphia, Pennsylvania, USA

Abstract Presented at The Diabetic Limb Salvage Meeting, March 2016

*The Effect and Feasibility of Low Frequency (22.5kHz) Contact Ultrasound Assisted Debridement (LFCUD) on Lower Extremity Wounds in a Vascular Surgery Clinic: The UltraHeal Pilot Study*

**Christine A. Murphy1,2 RN CETN(C) MCISc; Pamela Houghton1 PT PhD; Tim Brandys2 MD MEd FRCSC FACS; Gregory W. Rose2 MD MSc FRCPC**  
Western University1 and The Ottawa Hospital2, Canada

Ultrasonic Tangential Burn Excision Reduces Blood Loss  
**Abraham P. Houg MD MSE FACS, Sylvia J. Petrone MD FACS, Robin A. Lee MD, Christina Lee MD, Michael A. Marano MD FACS**  
Department of Surgery, Division of Burn Surgery Saint Barnabas Medical Center, Livingston, New Jersey  
Abstract Presented at the American Burn Association Meeting, April 2015

*Ultrasonic Debridement: An Effective and Superior Method of Wound Debridement*

**John M. Hiebert MD, FACS. FACN – Director Wound Care Centers, Saint Luke's Health System, Clinical Professor Plastic Surgery, University of Kansas**  
Abstract presented at Spring SAWC 2013

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© 2021 Biomedica Italia S.r.l. Diritti riservati
Biomedica Italia S.r.l. - Società unipersonale