WOUND CARE

SonicOne[®] O.R.

Debrider e Osteotomo ad Ultrasuoni

SonicOne® O.R.

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1. "A Closer Look At Ultrasonic Debridement" Issue Number: Volume 23 - Issue 8 - August 2010 (/issue/2235), Author(s); Martin E. Wendelken, D.P.M., RN, Lee Markowitz, D.P.M., and Oscar M. Alvarez, PhD"
2. "Split Thickness Skin Graft Incorporation after Ultrasonic Debridement for Wound Bed Preparation in Diabetic Foot Wounds." Sean M. Betesh, D.P.M., Michael I. Gazes, DPM, MPH, Peter Blume, DPM, FACFAS.
3. "Ultrasonic Tangential Burn Excision Reduces Blood Loss." Abraham P. Houg M.D., Sylvia J. Petrone M.D., Robin A. Lee M.D., Christina Lee M.D., Michael A. Marano M.D.
4. "Surgical Debridement with SonicOne; an initial look at the bacterial load pre and post low frequency direct contact ultrasound." Daria Abolghasemi, D.O., Karen D. Szymanski, D.O., Michael Baruch M.D., Jamshed Zuberi M.D.
5. "Ultrasonic Debridement An Effective and Superior Method of Wound Management." John M. Hiebert M.D. Abstract presented at SAWC 2013.
6. "Ultrasonic Tangential Burn Excision Reduces Blood Loss." Abraham P. Houg M.D., Sylvia J. Petrone M.D., Robin A. Lee M.D., Christina Lee M.D., Michael A. Marano M.D.
7. "Applications of the Ultrasonic Bone Cutter in Spinal Surgery - Our Preliminary Experience." Al-Mahfoudh R, Qattan E, Ellenbogen JR, Wilby M, Barrett C, Pigott T. Br J Neurosurg. 2013, Jul 10.
8. "Can Ultrasound Debridement Facilitate Biofilm Removal from Diabetic Foot Ulcers?" Melinda Bowlby, D.P.M., Peter Blume, D.P.M. Podiatry Today. Volume 27 - Issue 8 - August 2014.
9. "Surgical Debridement with SonicOne; an initial look at the bacterial load pre and post low frequency direct contact ultrasound." Daria Abolghasemi, D.O., Karen D. Szymanski, D.O., Michael Baruch M.D., Jamshed Zuberi M.D.
10. "An overview of factors maximizing successful split-thickness skin grafting in diabetic wounds." Ryan J. Donegan, D.P.M., Brian M. Schmidt, D.P.M., Peter A. Blume. D.P.M. Diabetic Foot & Ankle.

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

STUDI CLINICI

Surgical Debridement with SonicOne® O.R.: an Initial Look at the Bacterial Load Pre- and Post- Low Frequency Direct Contact Ultrasound

Daria Abolghasemi, D.O.1, Karen D. Szymanski, D.O.2, Michael Baruch M.D.3, Jamshed Zuberi M.D.4

Affiliations

1. Resident Department of General Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey
2. Fellow Section of Plastic and Reconstructive Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey, USA
3. Chief Section of Plastic and Reconstructive Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey, USA
4. Trauma and Critical Care Surgeon, Department of Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey, USA
Abstract Presented at The Diabetic Limb Salvage Meeting, March 2016

Split Thickness Skin Graft Incorporation after Ultrasonic Debridement for Wound Bed Preparation in Diabetic Foot Wounds
Sean M. Betesh, DPM1, Michael I. Gazes, DPM, MPH1, Peter Blume, DPM, FACFAS2

1. Resident, Department of Podiatric Surgery, Yale New Haven Hospital, New Haven, CT, USA
2. Assistant Clinical Professor of Surgery, Anesthesia, and Orthopedics and Rehabilitation, Yale School of Medicine, New Haven, CT, USA

Abstract Presented at The Diabetic Limb Salvage Meeting March 2016

The Clinical Implications of a New Wound Debridement Device That Combines Low Frequency Direct Contact Ultrasound and Vacuum Aspiration

Mark S. Granick M.D.1, Michael Baruch M.D.2, Wayne J. Caputo M.D.3, Paul M. Glat M.D.4

Affiliations

1. Division of Plastic Surgery, Rutgers New Jersey Medical School, Newark, New Jersey, USA
2. Division of Plastic and Reconstructive Surgery, St Joseph's Hospital and Medical Center, Paterson, New Jersey, USA
3. Clara Maass Medical Center, Barnabas Health, Belleville, New Jersey, USA
4. Division of Plastic Surgery, Division of Burn Unit, St. Christopher's Hospital for Children, Philadelphia, Pennsylvania, USA

Abstract Presented at The Diabetic Limb Salvage Meeting, March 2016

The Effect and Feasibility of Low Frequency (22.5kHz) Contact Ultrasound Assisted Debridement (LFCUD) on Lower Extremity Wounds in a Vascular Surgery Clinic: The UltraHeal Pilot Study

Christine A. Murphy1,2 RN CETN(C) MCISc; Pamela Houghton1 PT PhD; Tim Brandys2 MD MEd FRCSC FACS; Gregory W. Rose2 MD MSc FRCPC
Western University1 and The Ottawa Hospital2, Canada

Ultrasonic Tangential Burn Excision Reduces Blood Loss
Abraham P. Houg MD MSE FACS, Sylvia J. Petrone MD FACS, Robin A. Lee MD, Christina Lee MD, Michael A. Marano MD FACS
Department of Surgery, Division of Burn Surgery Saint Barnabas Medical Center, Livingston, New Jersey Abstract Presented at the American Burn Association Meeting, April 2015

Ultrasonic Debridement: An Effective and Superior Method of Wound Debridement

John M. Hiebert MD, FACS. FACN – Director Wound Care Centers, Saint Luke's Health System, Clinical Professor Plastic Surgery, University of Kansas
Abstract presented at Spring SAWC 2013

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© 2022 Biomedica Italia S.r.l. Diritti riservati
Biomedica Italia S.r.l.
Società unipersonale - Soggetta ad attività di direzione e coordinamento di Addlife Development AB