





LORDOTIC AND ANATOMICAL SHAPE

Optimal fit to the vertebral end plates



IMPLANT MADE OF BLACKARMOR®

Reduced artifacts in postoperative imaging procedures





RELIABLE OSSEOINTEGRATION

Rough titanium coating

inspired by nature – built by icotec

ICOTEC CERVICAL CAGE

Factsheet




EXCELLENT IMAGING QUALITY

focused on bone – not metal

ICOTEC CERVICAL CAGE


icotec Cervical Cages – Lordotic


| Description | Depth | Width | Height | Reference number |
|--------------------------------|-------|-------|--------|------------------|
| icotec Cervical Cage, lordotic | 13 mm | 13 mm | 5 mm | 12-40-0513 |
| | | | 6 mm | 12-40-0613 |
| | | | 7 mm | 12-40-0713 |
| icotec Cervical Cage, lordotic | 13 mm | 16 mm | 5 mm | 12-40-0516 |
| | | | 6 mm | 12-40-0616 |
| | | | 7 mm | 12-40-0716 |
| | | | 8 mm | 12-40-0816 |



icotec Cervical Cages – Anatomic

| Description | Depth | Width | Height | Reference number |
|--------------------------------|-------|-------|--------|------------------|
| icotec Cervical Cage, anatomic | 13 mm | 13 mm | 5 mm | 12-41-0513 |
| | | | 6 mm | 12-41-0613 |
| | | | 7 mm | 12-41-0713 |
| icotec Cervical Cage, anatomic | 13 mm | 16 mm | 5 mm | 12-41-0516 |
| | | | 6 mm | 12-41-0616 |
| | | | 7 mm | 12-41-0716 |
| | | | 8 mm | 12-41-0816 |
| icotec Cervical Cage, anatomic | 13 mm | 18 mm | 5 mm | 12-41-0518 |
| | | | 6 mm | 12-41-0618 |
| | | | 7 mm | 12-41-0718 |
| | | | 8 mm | 12-41-0818 |



Information needed to use the device is made available in electronic form; current and previous versions can be downloaded in electronic form at ifu.icotec-medical.com (code = ) or can be requested by email or phone from icotec. On request, icotec will provide a paper version within seven calendar days at no charge. The electronic versions can be viewed with a freely available PDF reader (e.g., Adobe Acrobat Reader, which can be downloaded at www.adobe.com).



ifu.icotec-medical.com

ico | cœt

 icotec ag

Industriestrasse 12
9450 Altstätten
Switzerland
Phone: +41 71 757 00 00
Fax: +41 71 757 00 01
info@icotec.ch
www.icotec-medical.com

CE 0297

Not available for sale within the United States.

0130EN.2019-01